

FOR IMMEDIATE RELEASE

**NEWS**

**FCI OPHTHALMICS**

P.O. BOX 465

Marshfield Hills, MA 02051

Tel: 800-932-4202 or 781-826-9060


Fax: 781-826-9062

[info@fci-ophthalmics.com](mailto:info@fci-ophthalmics.com)


[www.fci-ophthalmics.com](http://www.fci-ophthalmics.com)

Contact: Anne Bohsack

March, 2006



Reusable Flexible Iris Retractors



Disposable Flexible Iris Retractors

FCI Ophthalmics introduces the Reusable & Disposable Iris Retractors to supplement the company's line of cataract products including the Morcher<sup>®</sup> Capsular Tension Rings, Cionni Capsular Tension Rings, Pupil Dilator and Mackool Capsular Support System.

The Flexible Iris Retractors are used to gently stretch and stabilize the iris for maximum visibility during surgery. Once the iris is hooked and retracted, these bright blue retractors are held securely in position by sliding the tiny adjustable silicone "donut" up to the limbus.

The Flexible Iris Retractors are made of blue polypropylene fiber and are easily inserted through 0.5mm incision making them ideal for temporary iris dilation.

The reusable iris retractors are packaged sterile, box of 6 in a steam autoclavable container. The disposable iris retractors are packaged sterile, box of 5.

For more information on the Reusable and Disposable Flexible Iris Retractors and all our cataract products, contact FCI Ophthalmics.