

FOR IMMEDIATE RELEASE

NEWS

FCI OPHTHALMICS

P.O. BOX 465

Marshfield Hills, MA 02051

Tel: 800-932-4202 or 781-826-9060

Fax: 781-826-9062

info@fci-ophthalmics.com

www.fci-ophthalmics.com

Contact: Anne Bohsack

October, 2006


FCI Ophthalmics introduces Morcher Ring Calipers

Newly approved to US market, the Morcher Ring Caliper is a sterile temporary surgical device that enables the surgeon to create a more precise anterior capsulorhexis.

The use of the Ring Caliper is appropriate for primary or secondary IOL exchange procedures and is recommended for any surgical procedure where a calibrated anterior capsulorhexis of a precise diameter is desired. A rhexis of about 4.75 mm can be achieved with a Ring Caliper Type 5 (inner diameter 5.0 mm) or a rhexis of 6.0 mm with the Ring Caliper Type 6 (inner diameter 6.0 mm).

Because the Ring Caliper is placed directly on top of the anterior capsule, the patient's eye movement does not interfere with the capsulorhexis. The Ring Caliper is easily inserted through the incision site and can be positioned with an Osher Manipulator or Kuglen Hook.

The Morcher Ring Caliper is particularly helpful with new technology IOLs which require more precise capsulorhexis. It is also a helpful tool when teaching cataract surgery to residents as they learn to make a capsulorhexis.



For more information contact:

Anne Bohsack @ 800-932-4202 or info@fci-ophthalmics.com

#