

FOR IMMEDIATE RELEASE



FCI

**FCI OPHTHALMICS**

P.O. BOX 465

Marshfield Hills, MA 02051

Tel: 800-932-4202 or 781-826-9060


Fax: 781-826-9062

info@fci-ophthalmics.com

www.fci-ophthalmics.com

Contact: Anne Bohsack

May, 2003



## DISPOSABLE VITRECTOMY KITS

FCI Ophthalmics has received the FDA approval to market in the United States its DISPOSABLE VITRECTOMY KITS for use in vitreo-retinal surgery. The kit includes a lens, a self-blocking sutureless infusion cannula, and a pre-calibrated knife.

The lens, available in planoconcave or 30° prismatic styles, is made of clear silicone for superior optical visualization and requires no suturing. It may be easily moved by grasping its outer lip with forceps.

The pre-calibrated knife assures a tight seal at the sclerotomy sites for instruments, cutters and the self-blocking sutureless infusion cannula.

The kit includes the planoconcave lens, the infusion cannula and the pre-calibrated knife. Available separately: infusion cannula with pre-calibrated knife or the lenses. All packaged 5 per box.

Also from FCI: a wide selection of scleral buckles for retinal detachment surgery